

## Certificate

PROFIBUS Nutzerorganisation e.V. grants to

### SIEMENS AG, DF MC

Frauenauracher Str. 80, 91056 Erlangen, Germany

the Certificate No: **Z10834** for the PROFINET IO Device:

Model Name: SINAMICS G120D CU240D-2 PN-F FO  
Revision: SW/FW: V4.7; HW: 2  
Identnumber: 0x002A; 0x0512  
GSD: GSDML-V2.3-Siemens-Sinamics\_G120D-20140212.xml  
DAP: IDD\_CU240D-FO-V4.7, 0x00020609

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/>	PNIO_Version	V2.2
<input checked="" type="checkbox"/>	Conformance Class	A, B, C
<input checked="" type="checkbox"/>	PNIO_Tester_Version	V2.2.4
<input checked="" type="checkbox"/>	Tester	an testlabs, Prague, Czech Republic, PN000-122, IRT000-122

This certificate is granted according to the document:

“Framework for testing and certification of PROFIBUS and PROFINET products”.

For all products that are placed in circulation by **July 17, 2020** the certificate is valid for life.

Karlsruhe, October 10, 2017

Board of PROFIBUS Nutzerorganisation e. V.



(Official in Charge)



(Karsten Schneider)



(Dr. Jörg Hähnliche)