

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

SIEMENS AG, SI EP R&D CP PMC 2

Werner-von-Siemens-Str. 48, 92220 Amberg, Germany

the Certificate No: **Z10505** for the PROFINET IO Device:

Model Name: SIMOCODE pro V PN
 Revision: SW/FW: V 2.1; HW: 14
 Identnumber: 0x002A; 0x0904
 GSD: GSDML-V2.3-SIEMENS-SIMOCODEproVPN-20180830.xml
 DAP: DAP1: SIMOCODE pro V PN, 0x00000010

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/> PNIO_Version	V2.2
<input checked="" type="checkbox"/> Conformance Class	C
<input checked="" type="checkbox"/> Optional Features	Advanced, IRT
<input checked="" type="checkbox"/> PNIO_Tester_Version	V2.2.4
<input checked="" type="checkbox"/> Tester	SIEMENS AG, Fürth, Germany PN195-2, IRT071-1

This certificate is granted according to the document:

“Framework for testing and certification of PROFIBUS and PROFINET products”.

For all products that are placed in circulation by **July 31, 2023** the certificate is valid for life.

Karlsruhe, April 06, 2020



(Official in Charge)

Board of PROFIBUS Nutzerorganisation e. V.



(Karsten Schneider)




(Dr. Jörg Hähnliche)