

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

SIEMENS AG, PD PA CI R&D 2

Östliche Rheinbrückenstr. 50, 76187 Karlsruhe, Germany

the Certificate No: **Z11257** for the PROFINET IO Device:

Model Name: SIMATIC MV420
 Revision: SW/FW: V7.0.2; HW: 7
 Identnumber: 0x002A; 0x0B08
 GSD: GSDML-V2.31-Siemens-MV420-20150925.xml
 DAP: DAP6_MV420_PIB: Ident profile V7.0, 0x00000600

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/>	PNIO_Version	V2.2
<input checked="" type="checkbox"/>	Conformance Class	B
<input checked="" type="checkbox"/>	PNIO_Tester_Version	V2.2.4
<input checked="" type="checkbox"/>	Tester	SIEMENS AG, Fürth, Germany PN375-1

This certificate is granted according to the document:

“Framework for testing and certification of PROFIBUS and PROFINET products”.


For all products that are placed in circulation by **October 29, 2018** the certificate is valid for life.

Karlsruhe, February 07, 2017


Board of PROFIBUS Nutzerorganisation e. V.



(Official in Charge)



(Karsten Schneider)

(K.-P. Lindner)